


Modes of Operation:

- Reflection NSOM allows images of the intensity of light reflected from the sample to be obtained. The sample is illuminated via a probe aperture during the scanning and reflected light is directed by a mirror system through the objective of the inverted microscope to the PMT (photomultiplier tube).
- Luminescence NSOM allows images of local luminescence intensity to be obtained. The sample is illuminated during scanning via a probe aperture and transmitted light is directed to the PMT through an inverted microscope objective and notch filter.
- Transmission NSOM allows images of the intensity of the light transmitted through the transparent sample to be obtained. The sample is illuminated via a probe aperture during the scanning and transmitted light is collected by an inverted optical microscope and directed to PMT.


Different Laser Wavelength Options	
Laser Option	Wavelength
Trestles 20/50/100	700-950 nm
Mavericks-65	1240-1270 nm
Tamarak-Er	1540-1560 nm

Scan type by sample; by probe (optional)	
Measuring head scan range	80x80x3.5 μm ($\pm 10\%$)
Min scanning step	0.009 nm
X,Y sample positioning	5x5 mm
Positioning resolution	5 μm
Sample size	up to $\varnothing 100 \times 10 \text{mm}$
X,Y Closed-loop stage	
X,Y range	100x100 μm
Residual non-linearity	Better than 0.2%
Resolution	2 nm
Repeatability	30 nm (typically), less than 0.2% from the full range
Max. normal load	2 kg
Illumination laser	
Pulse length	<50 fs
Output power	>100 mW
Repetition rate	> 80 MHz
Wavelength range	740-950 nm




NSOM lithography on positive photoresist made with 488nm Ar laser. Time of exposure for different spots 0.1-0.5s. The result of the lithography is measured in Shear Force Mode.

Scan size: 3.7x3.7 μm .




Linearized DNA plasmid on mica. Shear Force image mode.

Scan size 1250x1250x1.2nm.



Gold particles deposited on glass. Transmission NSOM image.

Scan size: 27x27 μm .



GaAs quantum dots. Reflection NSOM image made with 442 He-Cd laser.

Scan size: 7x7 μm .

DEL MAR PHOTONICS

4119 Twilight Ridge | San Diego, CA 92130 | Tel: (858) 876-3133 Fax: (858) 630-2376
E-mail: sales@dmphotonics.com URL: www.dmphotonics.com